

Scooters and Power Wheelchairs Accessible on Peterborough Transit



Power wheelchairs

MODEL	WIDTH	LENGTH	ACCESS?	TRANSPORTATION TIE DOWN BRACKETS AVAILABLE?
Invacare TDX Models	24" - 25.5"	34" - 35.25"	YES	YES
Invacare Storm Models	24.75" - 25"	28.5" - 31.5"	YES	YES
Invacare Power Tiger	22.25"	34"	YES	NO
Invacare Pronto	24" - 25.9"	34" - 39"	YES	YES M91 & M51 only
Invacare Nutron	20" - 26"	40.75"	YES	NO
Pride Jazzy	21.5" - 27.25"	27.75" - 46"	YES	YES with synergy seating only
Pride Jazzy 1170XL plus & 1650	30" - 30.375"	45.25" - 47"	NO	N/A
Pride Quantum	23" - 27.25"	33" - 46"	YES	YES with synergy seating only
Pride Quantum 1650	30.375"	47"	NO	N/A
Sunrise: Quickie	24.5" - 27.75"	39" - 43.5"	YES	YES
Maple Leaf: MP 1, 2 & 5	23.6" - 25.6"	41.3" - 42.9"	YES	YES
Shoprider: UL8W, 888WNLS, P424M & P424L	21.5" - 25"	30" - 41"	YES	YES
Golden Technologies: Alante & Compass Series	23" - 25.5"	32.5" - 40"	YES	NO
Nexidea: HS 2800 & HS 5600	24" - 24.8"	38.7" - 39"	YES	NO
Handicare: Alex	24" - 25.5"	32.5" - 34.5"	YES	NO
Permobil: Corpus Series	24" - 28.5"	32.5" - 43"	YES	YES
Permobil: Corpus X850	28"	49" - 57"	NO	NO
Rascal: 320, 326, 312 Turnabout, 314, 327, 710	21" - 27.5"	25.25" - 34.5"	YES	NO
Sunrise: Quickie Rhythm, Xplore, Xperience & Groove	22.5"	39" - 43.5"	YES	YES
Sunrise: Quickie X2, S-646, S-646SE, P-222 SE & Smart Seat	24.75" - 27.75"	39"	YES	YES x2 only



the home medical equipment specialists

329 McDonnell St, Peterborough

(705) 740-2650 • 1-800-667-0158

Scooters and Power Wheelchairs Accessible on Peterborough Transit



Scooters

MODEL	WIDTH	LENGTH	ACCESS?
Pride Revo, Go-Go, Maxima, 3WH Legend, Victory and Celebrity X	19.12" - 25.5"	36.5" - 47"	YES
Pride Wrangler, Hurricane, Legend XL, Celebrity XL & DX SP	24" - 31"	51.25 - 63.5"	NO
Invacare Lynx, Auriga and Leo	20" - 24"	39.8" 48"	YES
Invacare Meteor	25.5"	57.5"	NO
Maple Leaf MP3 and 4B	23.6" - 24"	44.1" 44.9"	YES
Maple Leaf MP4A and C	23.6" - 24.8"	52.8" - 56.7"	NO
Quingo Classic, Air and Plus	22" - 25"	42" - 48"	YES
Quingo Sport	25"	53"	NO
Shoprider Gk8, 9 & 10, 778S & 778EL, 888SL/SE & 888SLN 889SL Series	21.5" - 24.5"	38" - 49"	YES
Shoprider 889SBN, 778XLSN & 889XLSN	28" - 29"	56" - 63"	NO
Handicare Fortress 1700/2000	24" - 24.5"	44" - 49.5"	YES
Handicare Trophy and Winner	26" - 26.5"	55" - 55.1"	NO
Ranger Freedom, Solo, Safari and Spirit	21" - 26"	40" - 48.5"	YES
NexIdea: HS 320-570 Models	21" - 24"	43" - 49"	YES
NexIdea: HS 580 Model	24"	50"	NO
Golden Technologies: Buzzaround, Literider & Companion Series	21.25" - 24"	36.5" - 47.5"	YES
Golden Technologies: Avenger	24.6"	57.5"	NO
Rascal: Balance 600B, Liteway 214 & 224, Ultralite 370 & WeGo 250	21" - 25"	41" - 46"	YES
Rascal 655F	27"	54"	NO

NOTE:

- All power wheelchairs and scooters which are not accessible on Peterborough's regular transit system are accessible on Peterborough's Handivan.
- All models mentioned are current 2010 models, older models may have different demensions.
- All manual wheelchairs with a wheelbase less than 30" wide are accessible on Peterborough's regular transit system.
- All walkers and rollators are accessible on Peterborough's regular transit system.
- Most walkers and rollators today are capable of folding allowing for easier transport.
- Any powerchair with a wheel-base 30" wide or greater and any scooter over 51" long is not accessible on Peterborough's regular transit system.
- The maximum weight capacity of the regular bus and handivan bus ramps is 750lbs.



the home medical equipment specialists

329 McDonnell St, Peterborough

(705) 740-2650 • 1-800-667-0158